
**BegoPal+ is a palladium containing cobalt-chromium alloy for fabrication of porcelain fused to metal dental restorations. BegoPal+ is available as powder and is processed by selective laser melting.**

non-sterile

CE 0197

ISO 9693-1 / ISO 22674

Availability: REF 61550 = 1,0 g

**Standard analysis, % by weight (elements)**


Co	41	Pd	25
Cr	25	Mo	9

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## Directions for use

**Indications for use:** BegoPal+ is a porcelain-fused-to-metal (PFM) noble alloy. It is suitable for fabrication of crown and bridge restorations. BegoPal+ is available as powder and is processed by selective laser melting (SLM).

**For professional use only: Rx only**

**Contraindications:** No contraindications are known. However, unwanted biological reactions such as allergies to contents of the alloy or electrochemically based reactions may very rarely occur. In case of known incompatibilities and allergies to contents of the metallic material it should not be used.

**Warnings:** Metal dust is harmful to your health. When grinding and blasting use suitable air extraction system / ventilation at the workplace and breathing mask type FFP3-EN149:2001!

**Precautions:** This device is for single use only. Do not reuse the alloy. This may result in loss of proper mechanical properties. In case of occlusal or approximal contact of BegoPal+ with a different alloy electrochemically based reactions may very rarely occur. Safety and effectiveness in treatment of children or treatment of pregnant or nursing woman have not been established.

**Adverse reactions:** No adverse reactions are known. Nevertheless, the rare case of occurrence of individual reactions against single components of BegoPal+ can never be excluded completely. In this case, the application of BegoPal+ should not be continued.

**Prescription device:** Caution: Federal law restricts this device to sale by or on the order of a licensed dentist.

**Modelling:** Minimum metal thickness (after grinding): for ceramic veneering 0.3 mm, for acrylic veneering with retention pearls 0.3 mm. Shape groove in cervical and palatal area. Avoid sharp edges.

**Finishing:** For finishing use fine carbide tools.

**Ceramics:** Use ceramics according to DIN EN ISO 9693-1 with firing temperatures of up to approx. 980 °C / 1800 °F (e.g. Ceramco 3, Creation, IPS InLine, Noritake EX-3, Vintage, Vintage Halo). Also suitable: Ceramics with reduced firing temperature (e.g. VITA VMK Master) and normal thermal expansion. Always follow the ceramic manufacturer's instructions! Always blast the surface to be veneered with pure Korox® 250 (250 µm particle size / corresponding gauge number: 60 grit and 3 – 4 bar) and clean the frame thoroughly (steam clean or boil in aqua dest.). After cleaning, hold with artery clamps and refrain from touching.

**Oxide firing:** Not recommended. If oxide firing is carried out to check the surface (900 °C / 1650 °F / 5 minutes with vacuum), always blast the oxide again as stated above. Clean thoroughly (steam clean or boil in aqua dest.).

**Ceramic firing:** Always apply basic material in two firing operations. The first coating thin (washbrand), the second coating opaque. Wash off frame under running water before application of next ceramic coating. Allow to cool down normally after firing. Remove ceramics only mechanically. Hydrofluoric acid (HF) corrodes the metal frame.

**Acrylic veneering:** The corresponding instructions of the manufacturers must be followed when processing the veneering systems.

**Final work:** Blast visible metal surfaces with Korox® 50 (50 µm particle size / corresponding gauge number: 220 grit) e. g. inside crowns, and then blast-polish external surfaces with Perlblast®. This is of course dispensable for external surfaces when it has been completed before firing. Then continue rubber-polishing the external surfaces with BEGO rubber polishers, and after that finish-polish with BEGO-cobalt-chrome polishing paste (blue, REF 52310). Clean thoroughly (steam clean or boil in aqua dest.).

**Soldering:** Soldering prior to firing with the flame: Wirobond® solder (REF 52622) and Fluxsol flux (REF 52531). Soldering after firing in the furnace: WGL solder (REF 61079) and Minoxyd flux (REF 52530). Cool down normally.

**Laser welding:** Filler material: Wiroweld wire Ø 0.35 mm (REF 50003) or Ø 0.5 mm (REF 50005).

**Warranty:** Whether given verbally, in writing or by practical instructions, our recommendations for use are based upon our own experience and trials and can only be considered as standard values. Our products are subject to a constant continuous development. Therefore alterations in construction and composition are reserved.

**Limitation of Liability:** Except where prohibited by law, BEGO will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability. The device labeling meets the recommendations of FDA applicable guidance documents.

**SLM rapid prototyping devices that are compatible with BegoPal+:**

EOSINT M 270 and EOSINT M 280 (EOS GmbH, Krailling/Germany) as well as SLM 125HL and SLM 250HL (SLM Solutions GmbH, Lubeck/Germany).

Do not re-use	Manufacturer	Consult instructions for use	Caution

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